

METROHM AG CH-9101 Herisau (Switzerland)

pH Meter


744

Series 01 ...







8.744.1003

95.10 Ti/kg


Selection of measured quantity

	<p>Press <pH/mV/°C> key. Each keystroke changes the measured quantity.</p>
---	--


Meaning of the symbols in the display

Display	Meaning
	<p>"A key must be pressed", e.g. <ENTER> or input of a value and <ENTER>.</p>
	<p>"Immerse electrode in buffer solution for calibration." Flashing beaker means "change buffer".</p>
	<p>Display of the calibration parameter "slope", relative slope of the electrode.</p>
	<p>Display of the calibration parameter "pHas", asymmetry pH.</p>
	<p>Drift display. Measured value drifts towards higher or lower value.</p>
	<p>The displayed temperature is measured with a Pt 1000.</p>
<p>hold 9</p>	<p>Hold buffer contains 9 values. Flashing "hold 3" means "the hold buffer is being viewed; the third value is displayed".</p>

Input of numbers

 <p>Flashing dot</p>	<p>Press <>> key several times until the dot at the bottom right of the position you wish to change flashes. Set the desired number with the <^> key. Press <ENTER> to store the value.</p>
---	--

Exit


	<p>Exit from</p> <ul style="list-style-type: none"> - calibration (1-point calibration) - inquiries and entries - viewing the hold buffer - error messages
---	--



Configuration

Call-up: <2ND> <CONFIG>.
Advance inquiries with <CONFIG> or <ENTER>.

Display	Meaning	Input range
	Drift-controlled transfer of a value to the hold buffer. 1 = yes: Value is transferred when <HOLD> key has been pressed and the drift criterion is met. 0 = no. 0 leads to the next inquiry:	0, 1
	Transfer of a value to the hold buffer after elapse of a preset time interval. Time starts when <HOLD> is pressed. Transfer can be stopped and restarted with <HOLD>. 0: Immediate transfer when <HOLD> is pressed.	0 ... 199.9 min
	Beeper 0 = no 1 = yes	0, 1
	Display contrast	1 ... 4
	Time interval for automatic instrument shutdown after the last keystroke. 0 = no automatic shutdown.	0 ... 99 min
	Display of the program version	




Parameters

Call-up: <2ND> <PARAM>.
Advance inquiries with <PARAM> or <ENTER>.

Display	Meaning	Input range
	Temperature. If a temperature sensor is connected, the temperature is continuously updated.	-199.9 ... 199.9 °C
	Calibration parameter "slope", relative slope of the electrode.	0.1 ... 199.9 %
	Calibration parameter "pHas", asymmetry pH.	0.00 ... 14.00



+



Preselection of a buffer series

Call-up: With instrument switched off, press <pH CAL> and switch on instrument. Selection with <pH CAL>, then <ENTER>.

	<ul style="list-style-type: none"> 1 Metrohm buffers 2 DIN/NBS buffers 3 Fisher buffers 4 Merck/Riedel deHaën buffers 5 Ciba/Geigy buffers SP Special buffers 	The buffer series currently selected is marked, e.g. "- 1 -"

Instructions for Use for 744 pH Meter


Table of contents

1. Getting started	
1.1 Connecting electrodes	1
1.2 Measuring pH values	2
1.3 pH calibration	3
1.4 Measuring redox voltages and temperature	4
1.5 Adjusting the temperature	4
2. pH calibration	
2.1 Selecting a buffer series	5
2.2 pH calibration with special buffers	6
2.3 General information on pH calibration	7
3. Storing measured values, hold function	
3.1 Storing measured values	8
3.2 Viewing measured values	9
3.3 Deleting all measured values in the hold buffer	9
4. Preparations	
4.1 Switching on	10
4.2 Mounting in the wall holder	10
4.3 Setting the operating mode, configuration	11
5. Error messages, problems and troubleshooting	
5.1 Error messages	12
5.2 Measuring problems	13
5.3 Diagnosis	15
6. Appendix	
6.1 Technical specifications	21
6.2 Explanation of technical terms	23
6.3 pH values of buffer solutions	24
6.4 Warranty and certificates	26
6.5 Scope of delivery and ordering designations	29
6.6 Accessories	29
Index	33



1. Getting started

If you are starting operation with the 744 pH meter for the first time, it is advisable to read through section 4, pages 10 ff.

1.1 Connecting electrodes



1.2 Measuring pH values



- Check whether your electrode contains sufficient KCl solution (3 mol/L). It should be filled to a depth just below the fill hole. During measurement, the fill hole should remain open. Immerse the electrode in the analysis solution until the diaphragm is covered (app. 2 cm).
- Switch on the 744 pH Meter: press the <ON/OFF> key.
- Press the <pH/mV/°C> key once or more until "pH" appears in the display. The second line always shows the temperature as pH values are temperature dependent. If a Pt 1000 temperature sensor is connected, a thermometer symbol appears in front of the temperature which indicates that the displayed temperature has been measured. Without Pt 1000, you should enter the current measuring temperature if it is not already correct, see page 4. The triangle in the display remains lit up until the pH value is "stable".
- Wait for the triangle in the display to disappear, i.e. your measured value is stable.
- Take the electrode out of the analysis solution and rinse it thoroughly with dist. water. If need be, carefully dry the electrode with a paper towel.
- Perform the next measurement.

For good pH measurements, a pH calibration should be performed, see following page.

After the measurement, the glass electrode must be stored in KCl solution, $c = 3\text{ mol/L}$ (in the reference electrolyte). Close the fill hole with the cap.

Never store in water: The AgCl of the reference system dissolves better in solutions containing chloride than in pure water (the chloro complex is formed). If an electrode is stored in water, AgCl can precipitate in the diaphragm (and clog it)!


Never store dry: Otherwise, the gel layer of the glass membrane will dry out.

You will find details on pH measurements in Metrohm Application Bulletin 188. This is available free of charge from your Metrohm agency.

1.3 pH calibration

For pH measurements, a calibration is needed from time to time as the measuring properties of the electrodes change. Perform a calibration, e.g. daily before starting your measurements.

Two buffer solutions are required for this. If you are using Metrohm buffer solutions, you can immediately start the calibration. Other buffer solutions must first be preselected, see page 5.



- Immerse the electrode in the first buffer solution.

- Press the <pH CAL> key.

- The display shows the temperature.

The finger means that a key must be pressed.

To change the temperature, use the <»> and <^> keys. Select the digit you wish to change with the <»> key. You will recognise the selected digit by the flashing dot at the bottom right. Now change the number using the <^> key: Press the key repeatedly until the desired digit is displayed. When the correct value is set, press <ENTER>. If you have a Pt 1000 temperature sensor connected, the temperature does not concern you, it will be measured directly.

- Your first buffer solution is now measured.

After the measurement, the display prompts you to immerse the electrode in the second buffer solution: The buffer beaker flashes. Carry this out and press <ENTER>.

- Your second buffer solution is measured.

- After the calibration, the slope and asymmetry pH (pHas) of the electrode are briefly displayed and the instrument is ready for pH measurements.

If the pH meter comes to a stop on display of the slope or the asymmetry pH, the values are outside the limits, see page 7. If you then wish to accept the value, press <ENTER>. Otherwise, press <pH/mV/°C> to retain the old value.

You have now calibrated the electrode and you are ready to start pH measurements.


You can quit the calibration routine at any time by pressing the <pH/mV/°C> key.

The results of the calibration can be examined and changed with the key sequence <2ND> <PARAM>, see Short Instructions under "Parameters".

1.4 Measuring redox voltages and temperature

For redox voltage measurements you need a platinum or a gold electrode (see page 31) and for temperature measurements a Pt 1000 temperature sensor.

Redox voltages:




- Press the <pH/mV/°C> key once or more until "mV" appears in the display.

If you have a Pt 1000 temperature sensor connected, the temperature together with the thermometer symbol will be displayed in the second line.

The **triangle** in the display remains lit up until the measured value is "stable".

Temperature:



- Press the <pH/mV/°C> key once or more until "°C" appears in the display.


The thermometer symbol lights up to show that the temperature has been measured.

The **triangle** in the display remains lit up until the measured value is "stable".

If you measure pH values or redox voltages and have a Pt 1000 temperature sensor connected, the measured temperature always appears in the second line of the display together with the thermometer symbol.

1.5 Adjusting the temperature

The measuring temperature can be set with the key sequence <2ND> <PARAM> .



- Press <2ND> <PARAM> .

- The temperature can now be set, see page 3. Measuring and calibration temperatures are identical. If you set a temperature in the calibration procedure, this is also assumed to be the measuring temperature. Conversely, any temperature set here applies to the calibration routine.

- Confirm the set value with <ENTER> .

- Quit the inquiry with <pH/mV/°C> .


2. pH calibration

2.1 Selecting a buffer series

The instrument automatically recognises the buffers. As the pH values of different buffer series differ slightly, it is necessary to preselect the buffer series you wish to use in the subsequent calibration.

Metrohm buffers are preselected in the factory. If you wish to use a different buffer series, you have to select this.

Proceed as follows to select the buffer series:



- Switch the instrument off.
- Press and hold the <pH CAL> key while switching on the instrument (<ON/OFF> key).
- The display for buffer selection appears. The number specifies the buffer series and the dashes "-" indicate which series is currently selected.


The following buffer series can be selected:

- 1 Metrohm buffers
- 2 DIN/NBS buffers
- 3 Fisher buffers
- 4 Merck/Riedel deHaën buffers
- 5 Ciba/Geigy buffers

If you do not wish to use any of these buffers, select

SP Special Buffers. For the calibration procedure with special buffers, see page 6.


- Select the buffer series with the <pH CAL> key: Press the key repeatedly until the desired buffer series is displayed. Confirm the series with <ENTER>.



2.2 pH calibration with special buffers

When you have selected a buffer series, turn to page 3 for the calibration procedure.

The calibration procedure with special buffers (buffer preselection SP) differs slightly from the "normal" calibration routine as the pH values of the buffers have to be entered:



- Immerse the electrode in the first buffer solution.
- Press the <pH CAL> key.
- The display shows the temperature, which you enter (keys <»>, <^> and <ENTER>) or which has been measured.
- Enter the pH value of your first buffer solution (keys <»> and <^>) and press <ENTER>. Note that the pH values of your buffer solution are temperature dependent!
- When you have entered the value, the first solution is measured.
You are then prompted in the display to immerse the electrode in the second buffer solution. Carry this out and press the <ENTER> key.
- Enter the pH value of your second buffer solution and press <ENTER>.
- The second buffer solution is measured.
- On completion of the calibration, the slope and pHas (asymmetry pH) are briefly displayed.

2.3 General information on pH calibration

EXIT

1-point



slope 98.5 %
cal

6.94 pH
pHas
cal

CONFIG
2ND

PARAM
pH CAL



- The pH calibration remains stored even when the instrument is switched off.
- **Exit** from the calibration procedure is possible at any time with the <pH/mV/°C> key. If the <pH/mV/°C> key is pressed before measurement of the first buffer solution, the old calibration data are retained. On completion of the measurement of the first buffer solution, there is a 1-point calibration.
- In a **1-point calibration**, a new value of pHas is recorded, the old slope is retained.
- To **continue** the calibration procedure, the <ENTER> or <pH CAL> key can be used. Any new values entered will be accepted only if the <ENTER> key has been pressed.
- If the values for the **slope** and **asymmetry pH** (pHas) are outside the limits, the calibration procedure stops and the value in the display flashes. This is intended to draw your attention to the unusual values found in your pH calibration. If you nevertheless wish to accept the value, press <ENTER>. If you press <pH/mV/°C>, the new value will not be accepted and the old one remains stored.
Limits:
Slope: $90.0 \% \leq \text{slope} \leq 105.0 \%$
Asymmetry pH: $6.40 \leq \text{pHas} \leq 8.00$
- **Slope** and **pHas** (asymmetry pH) can be viewed and entered with the key sequence <2ND> <PARAM>, see Short Instructions under "Parameters".
- For good pH measurements, the two buffer solutions should be at the same **temperature**. You should also perform the calibration at or very close to the temperature at which you will subsequently measure the pH.

3. Storing measured values, hold function

Up to 9 value pairs – measured value and corresponding temperature – can be stored with the <HOLD> key.




3.1 Storing measured values

Measured values can be stored with the <HOLD> key. If the hold function is activated, the display shows, e.g. "hold 3", i.e. 3 pairs of measured values are stored in the hold buffer.


Depending on the setting, the measured values can be transferred to the buffer storage in several ways:

- **Immediately.**
- When the measured value is **stable**, i.e. the drift criterion is met (the triangle in the display disappears).
- **Periodically** at preset time intervals.

The type of measured value transfer is set with the key sequence <2ND> <CONFIG> :



- Press <2ND> <CONFIG> . The displays can be advanced with <ENTER> .
- Transfer **stable** measured values. "d" represents drift.
1 = yes, transfer stable measured values.
0 = no.
The value can be changed with the <^> key.
The entry of 0 leads to the next inquiry:
- Measured value transfer **periodically** at the preset time interval. The time interval is set with the <>> and <^> keys.
The <HOLD> key is used to start the process. The display shows "hold 0" until the set time interval has elapsed. After this, the first measured value pair is transferred, the display changes to "hold 1" and the time interval again allowed to elapse etc.
The measured value transfer can be interrupted with the <HOLD> key (display "hold") and then restarted with the same key.
- If both entries are set to 0, measured value transfer occurs **immediately** after the <HOLD> key has been pressed.
- Exit the entry with <pH/mV/°C> .




3.2 Viewing measured values

The measured values can be viewed with the <»> key from the last to the first and with the <^> key from the first to the last.
"hold X" flashes in the display.

Exit viewing of the measured values with <pH/mV/°C> .

3.3 Deleting all measured values in the hold buffer


All measured values in the buffer are deleted with the key sequence



4. Preparations

4.1 Switching on

The display "init" appears after switching on the pH Meter if the power supply was interrupted (power supply plugged out or main power switched off).




Display "init"

confirm with <ENTER>








4.2 Mounting in the wall holder

To have the 744 pH Meter ready to hand, it can be mounted in the 6.2051.020 Wall Holder, which can be fixed to the wall.



4.3 Setting the operating mode, configuration

The operating mode can be configured with the key sequence <2ND> <CONFIG> (see also Short instructions under "Configuration").



- Press <2ND> <CONFIG>. Press <ENTER> to move through the displays.
- Settings for the hold function, see page 8.
- Off/on switching of beep
Setting 0 or 1 with the <^> key.
1 = beep on
0 = beep off
- Setting the display contrast
Settings 1...4 with the <^> key.
1 = normal setting
4 = setting at high temperatures when the background becomes highly visible.
- Setting of the automatic shutdown time of the pH meter after the last keystroke.
Settings 0...99 min with the keys <>> <^>.
0 = no automatic shutdown.
- Display of the program version (display "P 10").

5. Error messages, problems and troubleshooting

5.1 Error messages

The pH-Meter switches off automatically if the hold buffer contains 9 values and if you press the <pH CAL> key. Delete the hold buffer with <2ND> <CLEAR HOLD>.

- E1 In calibration, the voltage of solutions of the second buffer differs from the measured value of the first buffer solution by less than 6 mV.
Exit: <pH/mV/°C>
Remedy: – Change buffer solution
 – Use fresh buffer solutions
 – Check electrode
- E2 After the error E1, the voltage difference is still less than 6 mV.
Exit: <pH/mV/°C>
Remedy: – Use fresh buffer solutions. Are **both** solutions OK?
 – Check electrode. Is the cable connection also OK?
- E3 Buffer could not be assigned.
Exit: <pH/mV/°C>
Remedy: – Preselect correct buffer series
 – Use fresh buffer solutions
 – Check electrode
- E4 The temperature difference between the two buffer solutions is larger than 2°C.
Exit: <pH/mV/°C>
- E5 The calibration temperature is lower than 0°C or higher than 99°C. The pH values for buffers are no longer defined outside these limits.
Exit: <pH/mV/°C>
- E6 You have a Pt 100 temperature sensor connected instead of a Pt 1000, Metrohm ordering number 6.1110.100; or use the combined pH glass electrode with built-in Pt 1000, Metrohm ordering number 6.0238.000.
Exit: <pH/mV/°C>

5.2 Measuring problems

The following table contains a list of the most frequent malfunctions. The corresponding maintenance steps for pH glass electrodes are listed on page 14. The maintenance steps for the diaphragm are also valid for redox electrodes.

Malfunction	Possible cause	Remedial
Display of the pH Meter drifts.	<ul style="list-style-type: none"> - Liquid junction potential on reference electrode not constant. - Loose contact. - Electrode not plugged in or defective cable. 	<ul style="list-style-type: none"> ▶ Use 3 M KCl solution or other suitable reference electrolyte. ▶ Rectify fault.
Display of the pH Meter is "sensitive to hand capacity".	<ul style="list-style-type: none"> - Reference electrode is not filled. - Reference electrode filled with water by mistake. - Diaphragm clogged. - Measurement in poorly conducting solution. 	<ul style="list-style-type: none"> ▶ Top up with 3M KCl or other electrolyte solution, free from air bubbles. ▶ Empty out water and fill with 3M KCl. ▶ Clean diaphragm. ▶ Add conducting salt or use electrode with sleeve diaphragm.
Sluggish establishment of pH.	<ul style="list-style-type: none"> - Dirty diaphragm. - Adsorption at glass membrane. 	<ul style="list-style-type: none"> ▶ Clean diaphragm. ▶ Service glass membrane.
Slope too low	<ul style="list-style-type: none"> - Diaphragm contaminated. - Adsorption at glass membrane. - Deswollen glass membrane after measurements in anhydrous solvents. - Old electrode. - Poor buffer solution. 	<ul style="list-style-type: none"> ▶ Clean diaphragm. ▶ Service glass membrane. ▶ Soak electrode in water between the measurements. ▶ Regenerate glass membrane. ▶ Use fresh buffer solutions.
Shows same value in pH 4 and 7 buffers.	<ul style="list-style-type: none"> - Buffer solutions ok? - Crack in the glass membrane. - Connector damp or dirty. 	<ul style="list-style-type: none"> ▶ Use fresh buffer solutions. ▶ Use new electrode. ▶ Dry or clean connector.

Care of pH glass electrodes

Cleaning the diaphragm

- After measurements in media with *low chloride concentration* (precipitated AgCl in the diaphragm → dark-coloured diaphragm): place electrode overnight in conc. NH_3 , rinse with water and renew reference electrolyte.
- After measurement in *sulphide-containing* media (Ag_2S in the diaphragm → dark-coloured diaphragm): Place electrode in freshly prepared, slightly acidic 7% thiourea solution. Then rinse with water and renew reference electrolyte.
- With *organic contaminants*: Place electrode for app. 5 minutes in chromosulphuric acid at 80°C , then rinse thoroughly and renew reference electrolyte.
- If the above measures do not help: carefully file down diaphragm with a diamond nail file. The outflowing electrolyte should be visible as a dark ring.

Care and regeneration of glass membrane

- After measurement in *non-aqueous* media: soak electrode in water between measurements.
- After measurement in *protein-containing* media: Immerse electrode for several hours in a solution of pepsin in hydrochloric acid (5% pepsin in $c(\text{HCl}) = 0.1 \text{ mol/l}$). Then soak thoroughly.
- Regenerating of glass membrane: Immerse glass membrane either for 1 min in a 10% solution of ammonium hydrogenfluoride (NH_4HF_2) or for a few seconds in 40% HF. After the etching rinse for app. 10 s in an HCl solution ($\text{H}_2\text{O}:\text{conc. HCl} = 1:1$). Rinse electrode with water and allow to stand for 24 h in the storage solution.
- Clean electrode with appropriate solvent after measurements in non-aqueous, contaminating media (e.g. lacquers).

5.3 Diagnostic instructions (for program P 10)

The 744 pH Meter is a very precise and reliable instrument. Thanks to its rugged construction it is virtually impossible for external mechanical or electrical influences to have an adverse effect on its functions.


Although the occasional fault in the instrument cannot be ruled out completely, malfunctions are much more likely to be due to problems with the electrodes (see page 13) or to incorrect operation or use.

It is thus advisable to localise any trouble with the aid of this quick and simple diagnostic routine. The user will then need to call METROHM service only if there is really a fault in the instrument. Also, he can give the service engineer much more precise information by referring to the numbered steps of the diagnostic program.

When asking for help, always state the series number on the bottom plate (S...., see p. 10) and program version (P XX, see p.11).

Procedure

- The steps must be carried out in sequence and compared with the reactions of the 744 pH Meter (shown indented). In the "yes" case, continue with the next instruction.
- If the instrument does not show the expected response ("no" case), the diagnostic step in question must be repeated so as to rule out any operating error. Consistently wrong responses, however, very probably indicate a defect.
- Steps marked with a double arrow (\Rightarrow) allow re-entry into the test routine in the event of repetitions, provided the display shows

 Otherwise repeat step 1 b).

- Pressing the <pH/mV/°C> key returns the instrument to the user program. To re-enter the diagnostic program, see preceding paragraph.
- Error indication: the display shows 'E' and a number.

Equipment needed

- Reference voltage generator, e.g. Metrohm pH Simulator 2.642.0010 (or any voltage source + accurate DVM, class 0.1 mV)
- Insulated, high impedance cable link, e.g. Metrohm 6.2108.060
- Pt 1000 simulator, or resistor switchbox, class 0.1 %, or resistor 1 k Ω /0.1 % and suitable short cable

➤ 1. a) **Prepare Instrument for diagnostic test**

Switch off pH Meter.

Detach all external connections (cables at rear).

b) Before switching on, press <2nd>, switch on pH Meter; hold down <2nd> key until 'd - 1' appears.

➤ 2. **Keypad test**

Press <2nd> (as often as needed) until appears,


d - 3

then press <hold/enter>

0.

Press the keys in the order of their associated key number (see Fig. 1):

The corresponding key number will be shown
(except with 'mode' → quit).



If a wrong key number is pressed, this is displayed automatically with 'Ex' (x = number of wrong key pressed).

Fig. 1 Keypad with key numbers

➤ 3. **Display test**

Press <2nd> (as often as needed) until appears,

d - 4

then press <hold/enter>

The display test runs automatically through the steps 1-10 shown in Fig. 2.

The test can be stopped at any time and started again with <hold/enter>.

At the end of the test, the display again shows 'd - 4'.





Fig. 2 Display test

➤ 4. **EEPROM test**

Press <2nd> (as often as needed) until appears,




then press <hold/enter>



The entire display 'd - 6' flashes several times while the test is running. If the test is successful, 'd - 6' reappears.

➤ 5. **Beep test**

Press <2nd> (as often as needed) until appears,



then press <hold/enter>

A short tune is heard.

➤ 6. **Check Instrument calibration**

<pH/mV/°C>

The display shows the last method selected before switching off.

6.1 Check potential measurement

Press <pH/mV/°C> (as often as needed) until the display shows the mV range.

Use a screened, high-impedance insulated cable to connect a voltage calibrator (mV generator, pH Simulator, etc.) of class 0.1 mV - or a less precise source connected in parallel with an accurate voltmeter (0.1 mV) - to measuring input 'pH/mV' of the 744. Set voltage to <2000 mV and compare with displayed reading (tolerance ± 1 mV).

Check high impedance of input:
(if the generator allows this). Switch source to 'high impedance' ($R_i \approx 1000 \text{ M}\Omega$) and compare display with the value read off previously. The value must not differ by more than 1 digit.

Disconnect calibrator.

6.2 Check temperature measurement

Press <pH/mV/°C> (as often as needed) until the display shows '°C'.

Attach Pt 1000 simulator or resistor switchbox (or $1 \text{ k}\Omega/0.1 \%$) to socket 'Pt 1000' and read off the temperature ($1 \text{ k}\Omega$ gives $0 \text{ }^\circ\text{C}$, tolerance $\pm 0.1 \text{ }^\circ\text{C}$).

Disconnect Pt 1000 simulator (or switchbox).


End of diagnosis

The diagnosis carried out so far, together with the tests that run automatically (battery test, RAM and ROM tests), have checked the functions of the 744 pH Meter. For more detailed investigations or more selective repeats one can also work through steps 7-10.

» 7. Power supply check

The 744 pH Meter is not fitted with batteries. Consequently, in this test only the voltage of the power supply is displayed.

Press <2nd> (as often as needed) until appears,



A parallelogram-shaped display showing the text "d - 2".

then press <hold/enter>



A parallelogram-shaped display showing the text "X.XX" on the top line and "bAt" on the bottom line.

Example




A parallelogram-shaped display showing the text "5.47" on the top line and "bAt" on the bottom line.

The supply voltage should not fall below 3.5 V. If there is an insufficient voltage: Check power supply and replace if necessary.


To exit, press <pH/mV/°C> ('d - 2' reappears)

➤ 8. **ROM test**


Press <2nd> (as often as needed) until appears,



then press <hold/enter>




The test runs automatically. If there are no defects, three '-' appear on the display from left to right.




➤ 9. **RAM-Test**


Press <2nd> (as often as needed) until appears,



then press <hold/enter>




The test runs automatically. If there are no defects, three '-' appear on the display from left to right.



➤ 10. **Initialise RAM**

It may happen on rare occasions that severe interference upsets the processor functions, causing the system to crash. After such an event the RAM area has to be initialised. Although the instrument's basic data are retained, RAM initialisation should be done only when necessary, as the stored user data (pH calibration data, selected buffers, configurations, etc.) are then erased.


Execute step 1.



Press <hold/enter>



If initialisation cannot be avoided, press <hold/enter>



Otherwise quit with <pH/mV/°C>.

During initialisation three '-' appear on the display from left to right. The standard parameters are now loaded.

The data lost from the user memory then have to be entered again, see Instructions for use, Appendix.

Execute step 11.

➤ **11. Restore the original situation**

Reconnect the electrode and temperature probe detached at the start of the diagnostic routine, and carry out a brief function test with these.

Summary of key assignments in diagnostic tests

(before switching on, press <2nd>)

For repeat observations and special applications it may be an advantage to perform a particular check directly. The numbers for each test are therefore given below.

Press <2nd> as often as required.

d - 1 RAM initialisation - only when necessary!

d - 2 Battery check (power supply)

d - 3 Keypad test

d - 4 Display test

d - 5 RAM test

d - 6 EEPROM test

d - 7 ROM test

d - 8 Beep test

6. Appendix

6.1 Technical specifications

Measuring range

pH value

pH = 0.00 ... 14.00

Voltage

U = -1990 ... +1990 mV

Temperature

t = -130.0... +199.9 °C

Resolution

pH value

Δ pH = 0.01

Voltage

Δ U = 1 mV

Temperature

Δ t = 0.1 °C

Errors (without measuring probes)

pH value

\pm 0.01

over operational voltage range

\pm 0.14% of (meas.value-pH7)

over temperature of use range

\pm 0.20% of (meas.value-pH7)

Voltage

\pm 1 mV

over operational voltage range

\pm 0.05% of meas.value

over temperature of use range

\pm 0.18% of meas.value

Temperature

\pm 0.1 °C

over operational voltage range

\pm 0.1 °C

over temperature of use range

\pm 0.4% of final value

Measuring rate

pH value with alternating

1.9 measurements/s each

temperature measurement

Voltage with alternating

1.9 measurements/s each

temperature measurement

Temperature

3.8 measurements/s

pH measurement

Temperature compensation of the slope
in the temperature range

0.0...100.0 °C

pH calibration

Type

1- or 2-point-calibration

Buffer solutions which can be used for
the automatic buffer recognition,
temperature dependence of the pH
values considered automatically

Metrohm, DIN/NBS, Fisher,
Merck/Riedel deHaën, Ciba/Geigy

Special buffers without automatic
recognition

Signal amplifier

Input resistance

$> 10^{13} \Omega$

Offset current

$< 3 \cdot 10^{-13} \text{ A}$

Display	
Type	LCD
Digit height	10 mm for measured value 6.9 mm for temperature
Dialog	with pictograms
Housing	
Materials	
Housing	Polybutylenterephthalate (PBTP)
Keypad membrane	Polyester (PET)
Protection against static discharges	inner walls of housing coated with aluminum
Ambient temperature, humidity	
Operational temperature	-10...55 °C
Storage, transport	-40...60 °C
Humidity	5...85 % rel. humidity
Reference conditions	
Adjusting interval	annually
Temperature	23 ± 2 °C
Operational voltage	6.0 ± 0.1 V
Error influence quantiles	
Operational voltage range	6 V ± 5%
Temperature of use range	0...55 °C
Safety specifications	Designed and tested in accordance to IEC publication 1010, safety class III. This manual contains some information and warnings which have to be followed by the user to ensure safe operation and to retain the apparatus in safe condition.
Power supply	
Power supply	6.2117.050 or 6.2119.050
Current consumption	6V DC ± 5%, 01A max. 15mA with Pt1000
Dimensions	
Width	172 mm
Length	188 mm
Height	40 mm
Weight	
incl. stand rod	900 g

6.2 Explanation of technical terms

Asymmetry pH

See pH calibration.

Buffer solutions

Buffer solutions have an exactly known pH value at a specified temperature. They are used for the pH calibration.

Calibration

See pH calibration

Diaphragm

Part of the electrode which forms a conductive connection between the analysis solution and the reference electrolyte solution.

pH electrodes have either a round ceramic diaphragm directly above the glass bulb or a sleeve diaphragm.

Drift

Measured values are usually not immediately *stable* when an electrode is immersed in the analysis solution. The measured value changes somewhat over a period of time. This change is called drift. In the 744 pH Meter, such drifting values are marked in the display with a triangle. The triangle fades as soon as the measured value is more stable than a preset drift threshold.

Drift threshold:

pH: 0.028 pH/min


U/mV: 1.875 mV/min

t/°C: 0.974 °C/min

pH_{as}

See pH calibration

pH calibration



In the pH calibration, the measured voltage is plotted against the specified pH value of the buffer solution. This results in a straight line (at least in the middle pH range).

The intersection point of this straight line with the pH axis gives the *asymmetry pH*, pH_{as} .

The *slope* of the line at 25°C is theoretically 59.16 mV per pH unit. This slope is specified as a relative slope 1 = 100%. Real pH calibrations usually result in a somewhat lower slope.

pH value	The pH value specifies the acidity or alkalinity of a solution. At pH = 7, the solution is neutral. Solutions with pH values less than 7 are acidic, those with pH values greater than 7 alkaline.
Redox measurements	Voltage measurement with a platinum or a gold electrode. The ratio of a reduced to an oxidised substance is measured.
Reference electrolyte	The reference electrolyte solution is an essential part of the reference electrode. In most cases, KCl solution, 3 mol/L is used.
Slope	See pH calibration

You will find additional information in the Metrohm monograph "Electrodes in potentiometry" and in Metrohm Application Bulletin 188. Both these publications can be obtained free of charge from your Metrohm agency.

6.3 pH values of buffer solutions

Metrohm buffers

Buffer t/°C	pH 4.00 ±0.02	pH 7.00 ±0.02	pH 9.00 ±0.02
0	3.99	7.11	9.23
5	3.99	7.08	9.18
10	3.99	7.06	9.13
15	3.99	7.04	9.08
20	3.99	7.02	9.04
25	4.00	7.00	9.00
30	4.00	6.99	8.96
35	4.01	6.98	8.93
38	4.02	6.98	8.91
40	4.02	6.98	8.90
45	4.03	6.97	8.87
50	4.04	6.97	8.84
55	4.06	6.97	8.81
60	4.07	6.97	8.79
65	4.09	6.98	8.76
70	4.11	6.98	8.74
75	4.13	6.99	8.73
80	4.15	7.00	8.71
85	4.18	7.00	8.70
90	4.20	7.01	8.68
95	4.23	7.02	8.67

pH 4: potassium hydrogen phthalate
pH 7: potassium/sodium hydrogen phosphate
pH 9: borax

DIN/NBS buffers

Buffer t/°C	C	D	F
0	4.01	6.98	9.46
5	4.00	6.95	9.40
10	4.00	6.92	9.33
15	4.00	6.90	9.28
20	4.00	6.88	9.23
25	4.01	6.87	9.18
30	4.01	6.85	9.14
35	4.02	6.84	9.10
38	4.03	6.84	9.08
40	4.03	6.84	9.07
45	4.04	6.83	9.04
50	4.06	6.83	9.01
55	4.07	6.83	8.99
60	4.09	6.84	8.96
65	4.11*	6.84*	8.94*
70	4.13	6.85	8.92
75	4.14*	6.85*	8.90*
80	4.16	6.86	8.89
85	4.18*	6.87*	8.87*
90	4.21	6.88	8.85
95	4.23	6.89	8.83

C: potassium hydrogen phthalate
D: phosphate
F: borax
according to DIN 19266 (1979)

Fisher buffers

Buffer t/°C	pH 4.00 rot	pH 7.00 gelb	pH 10.00 blau
0	4.01	7.13	10.34
5	3.99	7.10	10.26
10	4.00	7.07	10.19
15	3.99	7.05	10.12
20	4.00	7.02	10.06
25	4.00	7.00	10.00
30	4.01	6.99	9.94
35	4.02	6.98	9.90
38	4.02*	6.98*	9.87*
40	4.03	6.97	9.85
45	4.04*	6.97*	9.81*
50	4.06	6.97	9.78
55	4.07*	6.97*	9.74*
60	4.09	6.98	9.70
65	4.11*	6.99*	9.68*
70	4.13*	7.00*	9.65*
75	4.14*	7.02*	9.63*
80	4.16*	7.03*	9.62*
85	4.18*	7.06*	9.61*
90	4.21*	7.08*	9.60*
95	4.23*	7.11*	9.60*

pH 4: potassium hydrogen phthalate (SB 101)

pH 7: potassium dihydrogen phosphate/NaOH (SB 107)

pH 10: potassium borate/carbonate/KOH (SB

Merck/(Riedel-deHaën¹) buffers

Buffer t/°C	pH 4.00	pH 7.00	pH 9.00
0	4.05	7.13	9.24
5	4.04	7.07	9.16
10	4.02	7.05	9.11
15	4.01	7.02	9.05
20	4.00	7.00	9.00
25	4.01	6.98	8.95
30	4.01	6.98	8.91
35	4.01	6.96	8.88
38	4.01*	6.96*	8.86*
40	4.01	6.95	8.85
45	4.00*	6.95*	8.82*
50	4.00	6.95	8.79
55	4.00*	6.95*	8.76*
60	4.00	6.96	8.73
65	4.00*	6.96*	8.71*
70	4.00	6.96	8.70
75	4.00*	6.96*	8.68*
80	4.00	6.97	8.66
85	4.00*	6.98*	8.65*
90	4.00	7.00	8.64
95	4.00*	7.02*	8.63*

pH 4: sodium citrate/sodium chloride

pH 7: potassium/sodium dihydrogen phosphate

pH 9: boric acid/KCl/NaOH

1): The values indicated by Riedel deHaën may vary until up to Δ pH 0.02.

Ciba/Gelgy buffers

Buffer t/°C	pH 4.00 (P01)	pH 7.00 (P10)	pH 9.00 (P12)
0	4.01	7.11	9.20*
5	4.00	7.08	9.15
10	4.00	7.05	9.10
15	4.00	7.02	9.05
20	4.00	7.00	9.00
25	4.01	6.98	8.96
30	4.01	6.97	8.91
35	4.02	6.96	8.88
38	4.03	6.95*	8.85
40	4.03	6.95	8.84
45	4.04	6.94	8.80
50	4.06	6.94	8.77
55	4.07	6.93	8.74
60	4.09	6.93	8.71
65	4.11*	6.93*	8.69
70	4.13	6.94	8.67
75	4.14*	6.94*	8.65
80	4.16	6.95	8.63
85	4.18*	6.96*	8.61
90	4.21	6.97	8.60
95	4.23	6.98*	8.59

pH 4: potassium hydrogen phthalate

pH 7: potassium/sodium hydrogen phosphate

pH 9: borax/potassium dihydrogen phosphate

* : interpolated or extrapolated values

The values without * correspond to the manufacturer's specifications.

6.4 Warranty and certificates

Warranty

The warranty regarding our products is limited to rectification free of charge in our workshops of defects that can be proved to be due to material, design or manufacturing faults which appear within 12 months from the day of delivery. Transport costs are chargeable to the orderer.

For day and night operation, the warranty is valid for 6 months.

Glass breakage in the case of electrodes or other glass parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts of outside manufacture insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With regard to the guarantee of accuracy, the technical specifications in the Instructions for Use are authoritative.

Concerning defects in material, construction or design as well as the absence of guaranteed features, the orderer has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases METROHM from any liability to pay compensation.

If any instruments and parts have to be returned, the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material, the parts must be packed in a dustproof package (for instruments, use of a plastic bag is imperative). If open assemblies are enclosed in the scope of delivery that are sensitive to electromagnetic voltages (e.g. data interfaces etc.) these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with built-in voltage source belong in a non-conductive protective packaging). For damage which arises as a result of non-compliance with these instructions, no warranty responsibility whatsoever will be accepted by METROHM.

Certificate of Conformity and System Validation

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity:	744 pH Meter
System software:	Stored in ROMs
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland

This Metrohm instrument has been built and has undergone final type testing according to the standards:

<i>Electromagnetic compatibility: Emission</i> EN55022 / class B, EN50081-1/92	Generic emission
<i>Electromagnetic compatibility: Immunity</i> EN50082-1/92	Immunity
IEC801-2/91 class 3, EN50082-2, NAMUR	Static discharge
IEC801-3	Radiated rf electromagnetic field immunity
IEC801-4	EI. fast transient requirements
IEC801-5	"Surges" immunity
EN50093	Voltage dips, short interruptions
<i>Security specifications</i> IEC1010, EN61010 class 2, UL3101-1	


The technical specifications are documented in the instruction manual.

The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality, analytical performance, and accuracy of results.

The features of the system software are documented in the instruction manual.

Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.

Herisau, Sept. 14. 1995



Dr. J. Frank
Development Manager



Ch. Buchmann
Production and
Quality Assurance Manager

Ionenanalytik • Analyse des ions • Ion analysis • Análisis iónico
744 pH Meter



Metrohm Ltd.
CH-9101 Herisau
Switzerland
Phone + 41 71 53 8585
Fax + 41 71 53 89 01



EU Declaration of Conformity

The METROHM LTD. company, Herisau, Switzerland hereby certifies, that the instrument:
744 pH Meter
meets the requirements of EC Directives 89/336/EWG and 73/23/EWG.

Source of the specifications:

EN 50081-1	Electromagnetic compatibility, basic specification Emitted Interference
EN 50082-1	Electromagnetic compatibility, basic specification Interference Immunity
EN 61010	Safety requirements for electrical laboratory measurement and control equipment

Description of the instrument:

pH Meter for measurements of pH, voltages and temperature

Herisau, September 14, 1995

Dr. J. Frank
Head of engineering dept.

Ch. Buchmann
Head of production dept.
Management representative

6.5 Scope of delivery and ordering designations

744 pH Meter	2.744.0014
including the following accessories:	
1 Clamping ring 10 mm	6.2013.010
1 Stand rod 30	6.2016.050
1 Electrode Holder	6.2021.020
1 Power Supply, Euro 220 V...240V/6V DC	6.2117.050
1 Hexagon Screw	V.017.8016
1 Instructions for Use for 744 pH Meter	8.744.1003
744 pH Meter	2.744.0015
like 2.744.0014 but with another Power Supply:	
1 Power Supply, USA 100V...120V/6V DC	6.2119.050
744 pH Meter	2.744.0024
like 2.744.0014, additional with	
1 LL-Solitrode with built-in Pt1000	6.0228.000
744 pH Meter	2.744.0025
like 2.744.0015, additional with	
1 LL-Solitrode with built-in Pt1000	6.0228.000

6.6 Accessories, Electrodes

6.6.1 General

Pt 1000 temperature sensor	6.1110.100
Cable for Pt 1000 temperature sensor, length 1 m	6.2104.080
length 2 m	6.2104.110
Wall holder	6.2051.020

6.6.2 For pH measurements

Solutions:

Ready-to-use solutions, 500 ml

Buffer pH 4	6.2307.100
Buffer pH 7	6.2307.110
Buffer pH 9	6.2307.120

Sets buffer concentrate 50 ml each, total 250 ml buffer solution

Buffer pH 4 and pH 7 and 3M KCl solution	6.2302.010
Buffer pH 4, pH 7 and pH 9	6.2304.000

Which electrode do you use for your measuring problem?

Normally, the 6.0228.000 Solitrode or 6.0238.000 electrode, both with built-in Pt1000, can be used. If you have any problems with your pH measurements, the following list may help you finding the right electrode.

Fill LL electrodes with aqueous KCl solution only.

Together with the electrodes you need a separate 6.2104.020 electrode cable (1 m)

1): 6.0726.100 separate reference electrode: 6.2106.020 cable (1 m); inner and outer electrolyte LiCl sat. in Ethanol.

Water

Boiler feed water (outer electr. 0.1 M KCl)	6.0219.100
Drinking water	6.0219.100, 6.0239.100
Ground water	6.0219.100, 6.0239.100
Rain water (outer electr. 0.1 M KCl)	6.0219.100
Swimming pool water	6.0219.100, 6.0239.100
Waste water (change electrolyte frequently)	6.0219.100, 6.0239.100

Foodstuffs

Beer	6.0239.100
Cheese (cleaning with pepsine, s.page 14)	6.0234.100, 6.0236.100
Cheese (< 15 °C; cleaning with pepsine, s.page 14))	6.0214.100
Cream (cleaning with pepsine, s.page 14)	6.0219.100, 6.0239.100
Fruit juices	6.0239.100
Fruits	6.0234.100, 6.0236.100
Fruits (< 15 °C)	6.0214.100
Jam (> 50 °C)	6.0219.100
Ketchup	6.0239.100
Mayonnaise (cleaning with pepsine, s.page 14)	6.0239.100
Meat (cleaning with pepsine, s.page 14)	6.0234.100, 6.0236.100
Meat (< 15 °C; cleaning with pepsine, s.page 14)	6.0214.100
Milk (cleaning with pepsine, s.page 14)	6.0219.100, 6.0239.100
Molasses	6.0239.100, 6.0219.100
Mustard	6.0239.100
New wines	6.0222.100
Salad dressing	6.0239.100
Vegetable juices	6.0239.100
Vine	6.0222.100
Vinegar	6.0222.100
Yogurt (cleaning with pepsine, s.page 14)	6.0219.100, 6.0239.100

Cosmetics, Detergents

Creams (non-aqueous, sep. meas./ref.electrodes)	6.0133.100 + 6.0726.100
Decalcifying agent	6.0222.100
Detergents	6.0219.100, 6.0235.100, 6.0239.100
Lotions (partly aqueous)	6.0219.100, 6.0239.100, 6.0235.100
Mouth wash	6.0239.100, 6.0235.100
Pastes	6.0239.100
Shampoo	6.0235.100, 6.0239.100
Skin (surface)	6.0217.000

Soap	6.0239.100, 6.0235.100
Tooth paste	6.0239.100
Tooth wash	6.0239.100, 6.0235.100
Washing agent, liquid	6.0235.100, 6.0239.100

Pharmaceutics, agrochemistry

Fertilizers (change electrolyte frequently)	6.0219.100
Ointments (non-aqueous, sep. meas./ref.electrodes)	6.0133.100 + 6.0726.100
Soil samples (slurry)	6.0219.100
Urine	6.0220.100, 6.0222.100

Leather, paper, textiles

Bleaching bath	6.0220.100, 6.0222.100
Dye bath	6.0239.100, 6.0219.100
Paper (surface)	6.0217.000
Tanning bath	6.0222.100, 6.0220.100
Textiles (surface)	6.0217.000

Lacquers, paints

Adhesives, alkaline (rinse with solvent)	6.0222.100
Lacquers (rinse with solvent)	6.0239.100
Paints, alkaline	6.0222.100
Suspensions (aqueous)	6.0239.100
Suspensions (non-aqueous, sep. meas./ref.electrodes)	6.0133.100 + 6.0726.100

Metals, galvanics

Caustic bath (change electrolyte frequently)	6.0219.100
Galvanic baths (change electrolyte frequently)	6.0219.100

Photo, printing

Photo baths (outer electrolyte KNO ₃ sat.)	6.0219.100
---	------------

General

Emulsions	6.0239.100
Low-ionic media (outer electrolyte 0.1 M KCl)	6.0219.100
Small sample volumes	6.0234.100, 6.0236.100
Small sample volumes (< 15 °C)	6.0204.100, 6.0214.100
Surface, e.g. skin, paper	6.0217.000

6.6.3 For redox measurement

Solutions:

Redox standard solution + 250 mV, 250 ml	6.2306.020
3M KCl solution, 250 ml	6.2308.020

Electrodes:

A separate electrode cable forms part of the electrode:

Length 1 m	6.2104.020
Length 2 m	6.2104.030
Length 3 m	6.2104.040

Comb. platinum electrode	6.0415.100
Comb. gold electrode	6.0413.100

Index

Texts which appear in the display are *printed in italics*. Keys are marked with < >.

A		G	
Accessories	29	Glass electrode	30
Asymmetry pH	23	- problems	13
B		H	
<i>bE</i>	11	Hand	3
Beaker	3	<i>hold</i>	8
Beep	11	<HOLD>	8
<i>bu</i>	5	I	
<i>bu 2</i>	3	<i>init</i>	10
Buffers		Input of values	3
- explanation	24	Installation	10
- pH values	24f	ISO	27
- selecting series	5	M	
C		Merck buffers	
<i>c</i>	11	- pH values	25
<i>cal</i>	3, 5ff	- selection	5
Calibration	3, 5ff	Measuring	2, 4
Certificates	27	- problems	13
Ciba/Geigy buffers		- temperature	4
- pH values	25	Measured values	
- selection	5	- storing	8
<CONFIG>	11	- viewing	9
D		- transferring	8
<i>d</i>	8	Metrohm buffers	
Diagnosis	15	- pH values	24
Diaphragm	23	- selection	5
DIN/NBS buffers		O	
- pH values	24	Ordering numbers	28
- selection	5	P	
Drift	23	<i>P</i>	11
Drift controlled value transfer	8	<PARAM>	4
Display contrast	11	pH	
E		- as	23
Electrodes	29	- electrodes	30
- connecting	1	- calibration	3, 5ff
- care and maintenance	14	- measurement	2
- selection	30	- value	23
Entering values	3	<pH CAL>	3, 5ff
Error messages	12	<pH/mV/°C>	2, 4
F		Program version	11
Finger	3	Pt 1000	1, 29
Fisher buffers			
- pH values	25		
- selection	5		

R

Redox

- electrodes	31
- measurements	4, 24
Reference electrolyte	1, 24
Riedel deHaën buffers	
- pH values	25
- selection	5

S

Shutdown time	11
Scope of delivery	29
Setting the operating mode	11
Slope	24
Special buffers	5, 6
Stable measured values	23

T

<i>t</i>	11
Technical specifications	21
Temperature	4
Temperature sensor	1, 29
Thermometer	2
Troubleshooting	13ff

V

Values

- enter	3
- viewing	9
- storing	8
Version	11

W

Wall holder	10
Warranty	26