

Instrument: ACQ-SM#C13UPA368M

Sample Manager: Wash Syringe Static Leak Test

Weak Wash Syringe Results:

Test Passed

Test Pressure: 190 psi

Compressed Volume: 29.2 μ L

Decompressed Volume: 38.0 μ L

Start Pressure: 194 psi

End Pressure: 193 psi

Pressure Changed: -1 psi

Average Pressure Leak: 0 psi/sec

Strong Wash Syringe Results:

Test Passed

Test Pressure: 190 psi

Compressed Volume: 28.2 μ L

Decompressed Volume: 30.0 μ L

Start Pressure: 196 psi

End Pressure: 191 psi

Pressure Changed: -5 psi

Average Pressure Leak: 0 psi/sec