


W1,44. 19 11 . 30

Pressure Test

G1312A: Bin Pump: DE33214455



| | Action | Result |
|----------|---|--------|
| | Prepare test as described in help (Press i to view) | Done |
| J | Close purge valve and press continue | Done |
| | Wait until pump is on | Done |
| | Pressure ramp up | Done |
| | Wait until pressure is stable | Done |
| | Measure pressure drop within 2 min | Done |
| | Pressure drop < 2.0 bar/min, meas.: 0.70 | Passed |
| | Cleanup test | Done |
| | Finished - Open purge valve to release pressure | Passed |

ed by: $\frac{\mathcal{R}arJolegn}{01-04-2019}$ Date: $\frac{\mathcal{I}}{\mathcal{I}}$ Executed by:

Signature: